

# DotDiver2.0

Automated LINE / Dot Immunoblot Analyzer



CE

## Product Highlights

- Fully Automated Immunoblot Analyzer
- Processing and Evaluation of LINE / Dot Immunoblots
- Comprehensive Product Portfolio for Autoimmune Diagnostics

## DotDiver2.0 – Immunoblot Analyzer for Autoimmune Diagnostics

**DotDiver2.0** is a compact benchtop immunoblot analyzer for fully automated processing and evaluation of LINE and Dot immunoblots to support the diagnosis of autoimmune diseases.

Barcoded **DotDiver2.0** test strips and cartridges ensure the correct combination of reagents and the quality assured handling of a wide range of immunoblots for rheumatology and gastroenterology in laboratory routine.

Ready-to-use **DotDiver2.0** reagents in pre-filled disposable cartridges minimize laborious hands-on time after manual sample application.

The user-friendly **DotDiver2.0** software guarantees the reliable evaluation of processed immunoblot test strips.

The powerful easy-to-use **DotDiver2.0** with a small footprint is indispensable for diagnostic routine services.



## DotDiver2.0 – Specifications



### Throughput

Samples	1 - 24
Assays	1 - 24

### Reagents

Test Strips	Barcoded
Cartridges	Pre-filled, barcoded

### Dimensions

Size (L x W x H)	340 mm x 420 mm x 365 mm
Weight	15 kg

### Connections

Ethernet / LIS	1 x
USB	1 x
HDMI	1 x

### Power

Requirements	24 V DC, 5 A
--------------	--------------




## DotDiver2.0

### AUTOMATED IMMUNOBLOT ANALYZER

- Fully automated LINE / Dot immunoblot benchtop analyzer with small footprint
- User-friendly software
- Simultaneous performance of up to 24 different tests
- Ready-to-use reagents and test strips
- Automated barcode identification of test strips and cartridges
- Integrated drying of processed test strips
- Evaluation of processed test strips
- Export of results in pdf- or xls-file format or as paper print out
- LIS connectivity
- Low-maintenance, no liquid handling

#### Product Information

### DotDiver2.0



#### Contact

GA Generic Assays GmbH

Ludwig-Erhard-Ring 3  
15827 Blankenfelde-Mahlow OT Dahlewitz  
Germany

Phone +49 (0) 33708 9286 0

Fax +49 (0) 33708 4417 25

info@genericassays.com

www.genericassays.com

#### Order Information

**DotDiver2.0**

**REF 5075**

