

DotDiver2.0

Automated LINE / Dot Immunoblot Analyzer



Product Highlights

- Automated Immunoblot Analyzer
- Processing and Evaluation of LINE / Dot Immunoblots
- Comprehensive Product Portfolio for Autoimmune Diagnostics

YOUR RELIABLE PARTNER IN AUTOIMMUNE DIAGNOSTICS

20 Years of Experience, 150 Partners in more than 100 Countries

DotDiver2.0 – Immunoblot Analyzer for Autoimmune Diagnostics

DotDiver2.0 is a compact benchtop immunoblot analyzer for automated processing and evaluation of LINE and Dot immunoblots to support the diagnosis of autoimmune diseases.

Barcoded **DotDiver2.0** test strips and cartridges ensure the correct combination of reagents and the quality assured handling of a wide range of immunoblots for rheumatology and gastroenterology in laboratory routine.

Ready-to-use **DotDiver2.0** reagents in pre-filled disposable cartridges minimize laborious hands-on time after manual sample application.

The user-friendly **DotDiver2.0** software guarantees the reliable evaluation of processed immunoblot test strips.

The powerful easy-to-use **DotDiver2.0** with a small footprint is indispensable for diagnostic routine services.



DotDiver2.0 – Specifications



Throughput

Samples	1 - 24
Assays	1 - 24

Reagents

Test Strips	Barcoded
Cartridges	Pre-filled, barcoded

Dimensions

Size (L x W x H)	340 mm x 420 mm x 365 mm
Weight	15 kg

Connections

Ethernet / LIS	1 x
USB	1 x
HDMI	1 x

Power

Requirements	24 V DC, 5 A
--------------	--------------



DotDiver2.0

Automated Immunoblot Analyzer

- Automated LINE / Dot immunoblot benchtop analyzer with small footprint
- User-friendly software
- Simultaneous performance of up to 24 different tests
- Ready-to-use reagents and test strips
- Automated barcode identification of test strips and cartridges
- Integrated drying of processed test strips
- Evaluation of processed test strips
- Export of results in pdf- or xls-file format or as paper print out
- LIS connectivity
- Low-maintenance, no liquid handling

Product Information

DotDiver2.0



Contact

GA Generic Assays GmbH

Ludwig-Erhard-Ring 3
15827 Blankenfelde-Mahlow OT Dahlewitz
Germany

Phone +49 (0) 33708 9286 0

Fax +49 (0) 33708 4417 25

info@genericassays.com

www.genericassays.com

Order Information

DotDiver2.0

REF 5075

